

Saliva Drug Screen

Multi-Test


The Saliva Drug Screen Device is a screening test for the simultaneous, qualitative detection of 7 drugs of abuse - **amphetamine, methamphetamine, cocaine, opiates (morphine), marijuana, benzodiazepines and oxycodone.**

First Sign Saliva Drug Screen Device for AMP/MAMP/COC/OPI/THC/BZO/OXY and their metabolites is a rapid oral fluid test that can be performed without the use of an instrument. The test utilises monoclonal antibodies to selectively detect elevated levels of specific drugs in human oral saliva.

Each of the tests is read individually and independently of one another.



The process is simple.



No line in the T means positive

SBM
Smith BioMed
NEW PLYMOUTH

www.smithbiomed.com