

Further Information

To find out more about any of the Drug and Alcohol Test Kits or any Point of Care Medical Diagnostic Tests refer: www.smithbiomed.com or talk to your pharmacist.

Alcohol and drug addiction is a chronic, often relapsing brain disease that causes compulsive use, despite harmful consequences to the addicted individual and to those around him or her.

Physical warning signs of drug abuse:

- Bloodshot eyes, pupils larger than usual.
- Changes in appetite or sleep patterns.
- Sudden weight loss or weight gain.
- Deterioration of physical appearance, personal grooming habits.
- Unusual smells on breath, body, or clothing.
- Tremors, slurred speech, or impaired coordination.



SBM
Smith BioMed

e. info@smithbiomed.com
p. +64 6 757 3550
f. +64 6 758 1539

19 Dawson Street
PO Box 8108
New Plymouth 4342
NEW ZEALAND

Helpline: 0508 biomed (0508 246 633)

TAPS No: PP8302



At home/industry/business PROFESSIONAL DRUG SCREENING TEST PACKS



- Use in the privacy of your own home
- Simple to use
- 99% accurate
- Available in single, dual and multiple drug test packs
- Complying with AS/NZS 4308:2008
- Always read the instructions and use as directed

Single Test Packs

Single tests for the detection of marijuana are available. These are a rapid, one step test for the qualitative detection of THC in urine.

Reliable - 99% accurate
Fast - Results in 5 minutes
Simple - One easy step



2 Test and 5 Test packs available.

Dual Test Packs

Dual tests for the detection of marijuana and methamphetamine (P) are also available. These are for simultaneous detection of marijuana and methamphetamine (P) in urine.



Multi Drug Integrated Test Cup

The Multi-Drug Test Cup is an immunochromatographic assay for rapid qualitative detection of drugs of abuse and their principle metabolites in urine at specified cut off concentrations.

It is a one step assay, detecting simultaneously for the six most common drugs of abuse found in Australasia. Each of the tests is read individually and independently of one another. The process is simple.

Features

- **Complying with AS/NZS 4308:2008**
- Detection of multiple drugs of abuse, conveniently integrated into a single test cup.
- Latest wicking technology has improved the speed and readability of the results.
- Improved fastening ensures a secure specimen collection.


Drug Types Detected

Abbreviation	Target Analyte	Cut-Off (ng per ml)
Amp.	Amphetamine (Speed)	300
Bzo.	Benzodiazepines	300
Coc.	Cocaine	300
Met.	Methamphetamine (P)	300
Mop.	Morphine	300
THC	Marijuana	50

To ensure the validity of the urine sample, the cup features a temperature strip and appropriate specimen validity tests, to avoid attempts to mask the results.

In addition to the above, Smith BioMed also has tests for:

- Synthetic cannabis
- Saliva/oral fluid
- Breathalysers



Multi Test